Tait Confidential Tait Electronics Ltd



TECHNICAL NOTE TN-998-SR

TM8000 Calibration Application Release v2.73

6 April 2005

Applicability

This Technical Note details the changes in the latest TM8000 Calibration Application v2.73.00.08.

1. Calibration Application Changes

Calibration App v2.73 Changes

- ➤ The TM8000 Calibration Application now has the ability to read the radio's serial number. The radio's Calibration file now defaults to be saved under the title "TM8000 Cal − 19xxxxxxx" when saving to disk. Raised as TIMS 33265.
- ➤ The TOOLS > DOWNLOAD... link to the TM8000 Download Application has been removed from the Calibration Application. This appeared in v2.71, but it did not provide the help or error pop-ups that are provided when accessed from the TM8100 PC Application, where it is still available. Raised as FOCUS 19432.
- Resolved an issue found in testing where recalibrating the RF Power on a radio and programming this back would produce an 'incorrect band' error. Raised as TIMS 43458.

NOTE: If a TM8000 requires re-calibration it is important to temporarily remove any programmed audio Tap-In or Tap-out settings. Found in PROGRAMMABLE I/O > AUDIO.

Depending on configuration these may inhibit the radio's audio path and consequently fail calibration testing.

Date of Print: 06/04/05

Tait Confidential Tait Electronics Ltd

Compliance Issues None.

CSO's – Please inform all logistics, sales and technical staff,

and dealers of this release.

2. Issuing Authority

Name and Position Graham Brenchley

of Issuing Officer Technical Support Engineer

Confidentiality Confidential – This message or document contains proprietary

information intended only for the person(s) or organisation(s) to whom it is addressed. All Recipients are legally obliged to not disclose Tait technological or business information to any persons or organisations without the written permission of

Tait.

Distribution Level Associate.

Document History Original Release 6 April 2005 GCB